



BIOTECH SUPPORT GROUP

AlbuSorb™

Albumin Depletion From Serum or Plasma


- Removes 30 mg albumin/ml, >90%
- Affinity-type equivalence, virtually no cross-reactivity with other proteins
- Disposable, no column regeneration or cross-contamination
- Economical new surface technology, not based on blue-dye or immuno-affinity chromatography
- Mild binding conditions maintains tertiary structure and simple transfer to secondary analysis
- The flow-through (unbound) fractions retain their enzymatic and biological activity
- Removes albumin from most species including human, sheep, bovine, mouse, goat, rat, and calf.

Poly-electrolytes are polymers with repeating units of stationary charges. AlbuSorb™ comes from a class of solid-phase, or surface-based, elastomeric poly-electrolytic surfaces that bind proteins through an empirically derived chemistry combining elements of polymer composition, cross-linking architecture and charge properties. As with bio-polymers like DNA and Heparin, governing their reactivity is the spatial presentation of the electrostatic groups along a flexible polymer chain.

Unlike immuno-affinity, the surfaces utilized are disposable eliminating cycle to cycle variance and cross-contamination. AlbuSorb™ is supplied as a powder. Simply weigh, centrifuge and/or filter, and recover the albumin depleted serum in the supernatant.

Cancer Sera Before and After AlbuSorb™

- 1: Normal pooled serum control
- 2: Flow-through from normal serum
- 3: Breast cancer pooled serum control
- 4: Flow-through from breast cancer serum
- 5: Lung cancer pooled serum control
- 6: Flow-through from lung cancer serum
- 7: Pancreatic cancer pooled serum control
- 8: Flow-through from pancreatic cancer serum



Gel Image: SDS-PAGE non-reduced,
Criterion™ Tris.HCl (Bio-Rad) 4-15%

Note: All samples are from human female ages 40-60



BIOTECH SUPPORT GROUP

| Product | Size | # Serum Preps | Item No. |
|------------------|---------|-------------------------------|----------|
| AlbuSorb™ | 1 gram | 20, 25 μ l Serum Samples | A185-1 |
| AlbuSorb™ | 6 grams | 120, 25 μ l Serum Samples | A185-6 |

| Items | Item No | Item No | Reagent |
|-----------------------------------|------------------|------------------|----------|
| AlbuSorb™ | A185-1 (1 grams) | A185-6 (6 grams) | Supplied |
| Binding Buffer BB1, pH 7.5 | 30 ml | 180 ml | Supplied |



| Typical Performance | AlbuSorb™ | AlbuSorb™ PLUS |
|--|---------------------------------------|--|
| Serum Sample Volume | 25 μ l | 25 μ l |
| Albumin Removal | >90% | >85% |
| Immunoglobulin Removal | - | >85% |
| Recovered Protein Mass | 500-600 μ g (Albumin depleted) | 400-500 μ g (Albumin + Ig depleted) |
| LC-MS/MS unique proteins (single 3 hr gradient) | 350-400 | 350-400 |



BIOTECH SUPPORT GROUP

AlbuSorb™ Applications

| Sample Type | Disease | Analysis |
|----------------------|----------------------|--------------|
| Rat serum | Cancer | MALDI |
| Rat serum | Diabetes | Western Blot |
| Human Synovial fluid | Rheumatoid Arthritis | 2DE |
| Human Urine Exosomes | Diabetes | LC-MS/MS SRM |

PROTOCOL – Based on processing 25 µl Serum

For best results – the lysate should be clear and free of colloidal material. We recommend first filtering through a 0.45 µm syringe-type filter before beginning the prep.

1. Weigh out 50 mg of AlbuSorb™ powder in a spin-tube/microfuge tube.
2. Add 400 µl of **Binding Buffer BB1** to condition the AlbuSorb™ powder. Shake it manually/vortex for 3 min and then centrifuge for 2 minutes at 3000 rpm. Discard the supernatant.
3. Repeat step-2
4. As a requirement for albumin binding, add 250 µl of the **BB1 Buffer** and then add 25 µl of the serum to **Step 3**. Mix for 10 minutes on a rotating shaker.
5. Centrifuge for 4 minutes at 10,000 rpm, **supernatant contains serum proteins minus albumin. Note – when observing proteins on SDS-PAGE (4-15%), other proteins migrate to the same region as Albumin, and may not be fully resolved.**
6. Optionally the pellet (**mostly albumin**) can be eluted with 200 µl of **stripping buffer (0.2M Tris + 0.5M NaCl, pH 10 by mixing on a shaker for 10 min)** and centrifuge for 4 minutes at 10,000 rpm.

The protocol can be scaled up or down proportionally to adjust for different serum volumes. The surface amount can be adjusted to accommodate more or less albumin removal.



BIOTECH SUPPORT GROUP

References

Cerebrospinal Fluid

Gwenael Pottiez, Pawel Ciborowski. [Proteomic Profiling of Cerebrospinal Fluid Expression Profiling In Neuroscience](#) Neuromethods.2012;64:245-270

Synovial fluid

Happonen KE, Fürst CM, Saxne T et al. [PRELP protein inhibits the formation of the complement membrane attack complex](#).Journal of Biological Chemistry.2012;287(11):8092-100

Serum

Holmberg R, Refai E, Höög A.[Lowering apolipoprotein CIII delays onset of type 1 diabetes](#). Proceedings of the National Academy of Sciences.2011;108(26):10685-9.

Tang MX, Ogawa K, Asamoto M. [Effects of Nobiletin on PhIP-Induced Prostate and Colon Carcinogenesis in F344 Rats](#) Nutrition and Cancer.2011;63(2):227-33

Holmberg, Rebecka [Apolipoprotein CIII and Ljungan virus in diabetes](#) 2010. Doctoral Thesis

Lu Q, Zheng X, McIntosh T [Development of different analysis platforms with LC-MS for pharmacokinetic studies of protein drugs](#). Analytical Chemistry.2009;81(21):8715-23

Urine

Zubiri, Irene, et al. [Diabetic nephropathy induces changes in the proteome of human urinary exosomes as revealed by label-free comparative analysis](#). Journal of Proteomics (2013).

Patent

Berggren, Per Olaf, Yang, Shao-Nian. 2012. [Methods For Treating And/Or Limiting Development Of Diabetes](#).U.S. Patent 20120328630 Kind Code: A1, filed June 25, 2012, and issued December 27, 2012.

CONTACT US

We welcome your questions and comments regarding our products.

Call 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.
 Fax 732-274-2899
 Email sales@biotechsupportgroup.com
 Mail 1 Deer Park Drive, Suite M, Monmouth JCT, NJ 08852, USA