



BIOTECH SUPPORT GROUP

AlbuVoid™

Albumin Depletion Plus Low Abundance Serum Protein Enrichment from Serum, Plasma, Tissues And Culture Media

- Albumin voids in flow-through >95%, with <30 minute bind/wash/elute protocol
- Low abundance enrichment equivalent or better than hexa-peptides or immuno-affinity
- Disposable, cost-effective, no column regeneration or cross-contamination
- Mild elution maintains tertiary structure and simple transfer to secondary analysis
- Low abundance protein and enzyme enrichment, suitable for all analytical platforms (i.e., LC-MS, Immunoassay, functional)
- Species agnostic, validated on human, sheep, bovine, goat, fish, etc.
- No molecular weight or pI bias
- For LC-MS, optional seamless On-bead protocols (BASP™) workflows and unique proteolytic efficiencies
 - No in-gel digests, no solution digests, no C18 desalting, more consistent, reproducible results
 - Compatibility with quantitative label (i.e., iTRAQ) and label-free LC-MS methods

AlbuVoid™ is an albumin depletion reagent kit, however the beads do not bind albumin. It removes albumin from serum and plasma samples while enriching low abundance proteins on the beads. The **AlbuVoid™** protocol uses mild buffers; the protocol conditions are very gentle so that native enzyme and functional activity is retained in elution fractions.

AlbuVoid™ beads are derived from the **NRicher™** platform chemistry; a porous silica-bead library of individual imperfect fit polymeric ligands. The library was designed to facilitate weak binding of proteins, allowing for preferential displacement of the stronger bias binding proteins, at or above the sub-proteome saturation of the beads.


The lower abundance enriched sub-proteome that binds to **AlbuVoid™**, can be eluted off without significant carry-over of albumin. It is ideal for applications involving discovery and targeted proteomics, enzyme assays, immunoassay and microarrays, 1D & 2D gel electrophoresis and LC-MS.

In addition, all **NRicher™ beads, including AlbuVoid™** have been adapted to a protocol specifically designed for LC-MS applications whereby the low abundance proteome adsorbed to the beads can be Trypsin processed to its peptide constituents. This is called Bead Assisted Sample Prep or BASP™; the protocol is included as an optional digest method.



For targeted proteomics, the **NRicher™** knowledgebase of over 2000 serum proteins is downloadable, and can help select the best product/method(s) for particular protein(s). Go to: <https://www.biotechsupportgroup.com/category-s/335.htm>



BIO|TECH SUPPORT GROUP




Two-dimensional Gel Electrophoresis Analysis of AlbuVoid™ treated sheep serum. Samples were reduced, alkylated and protein normalized (176 µg total protein per gel). Second dimension was on 8-16% polyacrylamide gels. Gels 1 & 2 were stained with Coomassie. The circled regions indicate the albumin zone. Gel 1: Sheep serum sample. Gel 2: AlbuVoid™ Eluate. Gel 3: Same as 2 - AlbuVoid™ Eluate but restained with SilverQuest (Invitrogen) silver stain. The differences between the gels illustrate the efficiency of albumin removal. The vast majority of the other serum proteins are retained, with no intrinsic bias with respect to molecular weight or pI.



Typical Performance	AlbuVoid™	AlbuVoid™ PLUS
Serum Sample Volume	100 - 200 µl	25 - 50 µl
Albumin Removal	>85%	>90%
Immunoglobulin Removal	-	>80%
Recovered Protein Mass	2000-4000 µg (Albumin depleted)	500-1000 µg (Albumin + IgG depleted)
LC-MS/MS observable proteins	700-1000	800-1200

SELECTIVE BINDING WORKFLOW





BIOTECH SUPPORT GROUP

Items Required	5 Prep (part of Catalog # AVS-05, AlbuTrial Kit)	10 Prep	50 Prep	Item
AlbuVoid™ Beads	0.25 gram	0.5 gram	2.5 grams	Supplied
Binding Buffer AVBB (0.05M HEPES, pH 6.0)	4 ml	8 ml	40 ml	Supplied
Wash Buffer AVWB (0.05M HEPES, pH 7.0)	5 ml	10 ml	50 ml	Supplied
Elution Buffer AVEB (0.25M Tris + 0.5M NaCl, pH 10.0-10.5)	2 ml	4 ml	20 ml	Supplied
Spin-filter & tube assemblies	5	10	50	Supplied

Additional Spin-Filters (low protein binding, 0.45 µm filter element) can be purchased separately, please inquire.

If there are any questions about compatibility or substitution with other buffers, please contact us.

PROTOCOL – Based on processing 100-200 µl Serum

For best results – the lysate should be clear and free of colloidal material. We recommend first filtering through a 0.45 µm syringe-type filter before beginning the prep. Depending upon the quality of the sample, centrifugation times can be adjusted to increase g's or time, sufficient to process the sample through the beads.

Product	Size	Total samples processed	Item No.
AlbuVoid™	5 Preps	5 x 200 µl samples	AVK-05
AlbuVoid™	10 Preps	10 x 200 µl samples	AVK-10

Note: Please contact sales@biotechsupportgroup.com for prices in bulk quantities.

The protocol can be scaled up or down proportionally to adjust for different volumes. The bead amount can be adjusted to accommodate more or less Albumin removal.

Albumin and IgG Removal Products

AlbuSorb™ PLUS and **AlbuVoid™ PLUS** are products that deplete both Albumin and IgG. For more information, go to:

<https://www.biotechsupportgroup.com/AlbuVoid-PLUS-p/np-avk.htm>

<https://www.biotechsupportgroup.com/AlbuSorb-PLUS-p/apk285.htm>

1. BEAD CONDITIONING. Weigh out 50 mg of **AlbuVoid™** beads in a spin-filter. Add 250 µl of **Binding Buffer AVBB**. Vortex or mix thoroughly for 5 minutes at room temperature followed by centrifugation for 2 minutes at 1,000 g's. Discard the filtrate. Repeat step-1.



BIOTECH SUPPORT GROUP

2. **SAMPLE PROCESSING.** Add 200 μ l of **Binding Buffer AVBB** to beads followed by 200 μ l of the Serum to the beads. Vortex or mix thoroughly for 10 min and then centrifuge for 4 minutes at 5,000 g's.
3. Remove the filtrate as Flow-Through, containing the unbound Albumin.
4. To the beads, add 500 μ l of **Wash Buffer AVWB**. Vortex for 5 min and centrifuge for 4 minutes at 5,000 g's. Discard the **Wash** filtrate.
5. Repeat Wash Step-4. The **AlbuVoid™** beads contain the enriched Albumin-depleted sub-proteome.
6. **BEAD ELUTION.** To the beads, add 400 μ l of **Elution Buffer AVEB**. Vortex for 10 min and centrifuge for 4 minutes at 5,000 g's. Recover the filtrate as elution (albumin depleted proteins). The proteome eluate (0.25M Tris + 0.5M NaCl, pH 10.0-10.5) is ready for further functional, proteomic or LC-MS analysis.

As an option for LC-MS sample preparation, the bead assisted on-bead digestion protocol (BASP™) is provided below. The digest buffer is not supplied but the Wash Buffer (0.05M HEPES, pH 7.0), can be made up, or comparable buffers (0.02-0.5M, pH 6-7) can be used. Higher pH buffers are not recommended.

6. Using **Wash Buffer**, prepare to 10mM of DTT concentration, and add 200 μ l to the beads and vortex for 10 minutes and incubate for 30 minutes at 60C.
7. Cool the samples to RT, add suitable volume of Iodoacetamide to 20mM and incubate in the dark for 45 minutes
8. Centrifuge for 4 minutes at 5,000 g's, and discard filtrate. Rinse the bottoms of the spin-filter tubes with 500 μ l of 50% ACN, **Wash Buffer** twice, to remove any traces of the filtrate.
9. Add 16 μ g trypsin in 200 μ l **Wash Buffer** to the beads and keep at 37°C for a minimum 4 hours to maximum overnight. Overnight is recommended to start with. In select targeted circumstances, 2 hours may be sufficient.
10. Centrifuge for 4 minutes at 5,000 g's, and retain digested peptides filtrate.
11. To further extract remaining peptides, add 300 μ L 10% formic acid, vortex 10 min, centrifuge for 4 minutes at 5,000 g's., and combine this volume with volume from step 10.
12. Total is about 500 μ l. Prepare to desired final concentration. Store at -80°C until LC-MS/MS.



BIOTECH SUPPORT GROUP

Selection of AlbuVoid™ References:

Secretome from Cell Culture (BSA Removal)

Paulsen, Bruna, et al. "[Reproducible differentiation of pure ovarian support cells from clinical-grade hiPSCs as a novel infertility treatment.](#)" *bioRxiv* (2024): 2024-04. To gain insight into the potential mechanism of action of these cells during oocyte maturation beyond transcriptomics readouts, the investigators performed secretome proteomics to investigate proteins that overexpressed in OSC after 24 hours in vitro in comparison with OSCs prior to culture (0 hour) from cells derived, from both CG-hiPSC and RUO-hiPSC. For proteomic analysis, the article states, "The **conditioned media ... and then passed through albumin depletion columns (AlbuVoid) to eliminate HSA-derived albumin.**"

Jenull, Sabrina, et al. "The histone chaperone HIR maintains chromatin states to control nitrogen assimilation and fungal virulence." *Cell Reports* 36.3 (2021): <https://www.sciencedirect.com/science/article/pii/S2211124721008196#mmc1>. The authors report a pivotal role for the HIR histone chaperone complex in modulating virulence of the human fungal pathogen *Candida albicans*. The article states for "Cell-free supernatants from 16 hours YNB-BSA (0.025% BSA) cultures grown at 30°C were used for Mass-Spec analysis. Collected supernatants were lyophilized and dissolved in 400 µl of water for **AlbuVoid™** treatment for albumin depletion...**Albumin-free enriched secretory proteome** was eluted from beads".

Serum/Plasma

Jing, Lun, et al. "[PROTEOMIC ANALYSIS IDENTIFIED LBP AND CD14 AS KEY PROTEINS IN BLOOD/BIPHASIC CALCIUM PHOSPHATE MICROPARTICLE INTERACTIONS.](#)" *Acta Biomaterialia* (2021). Here, in a LC-MS/MS proteomic study, the article describes use of **HemoVoid™** and **AlbuVoid™** prior to LC-MS analysis, "...After **albumin depletion**, analysis of the significant deregulated proteins showed that **27 signaling pathways significantly changed in blood cells...**"

Poillet-Perez, Laura, et al. "[Autophagy maintains tumour growth through circulating arginine.](#)" *Nature* (2018): 1.

One *in vivo* model to study autophagy is whole-body deletion of the essential autophagy gene *Atg7* in adult mice which causes a systemic metabolic defect that manifests as starvation intolerance. In order to measure the systemic proteomic response of such deletion in this study, **AlbuVoid™ was chosen as one method to remove albumin from mice and enrich the low abundance proteomes from serum.**

Vialaret, Jerome & Kadi, Sarah & Tiers, Laurent & O Flynn, Robin & Lehmann, Sylvain & Hirtz, Christophe. (2018). Albumin depletion of human serum to improve quantitative clinical proteomics. *Current Topics in Peptide & Protein Research* 19. 53-62. <http://www.researchtrends.net/tia/abstract.asp?in=0&vn=19&tid=26&aid=6192&pub=2018&type=3>

In this work, the investigators focused on depleting albumin from human serum samples using an albumin depletion and low abundance protein enrichment kit – AlbuVoid™, which enabled the detection of several low-abundance proteins. By employing an optimized protocol, enriched proteins known as biomarkers for various diseases were identified. The authors concluded that the **AlbuVoid™ depletion method proved to be faster and more cost-effective than antibody based methods.**

David L. Wang, Chuanguang Xiao, Guofeng Fu, Xing Wang and Liang Li. "[Identification of potential serum biomarkers for breast cancer using a functional proteomics technology.](#)" *Biomarker Research* (2017) 5:11.

The article states "The **most dramatic difference for enzyme activity detection in using the AlbuVoid™ for serum protein enrichment was demonstrated** in the case of protease activity analysis. Compared with the direct serum proteinase measurement, both the levels and species of proteases were increased significantly in the enriched serum sample. ..., **and it is necessary to use AlbuVoid™ to enrich these low level proteases to bring them to a high enough level to be detected.**"

For a full list of Albumin Removal references, visit:

<https://www.biotechsupportgroup.com/References-s/138.htm#albumin-depletion>



BIOTECH SUPPORT GROUP

For Targeted Proteomics

NRicher™ Bead Platform Provides Unique Sub-Proteome Biases And Fit For Purpose Opportunities for Targeted LC-MS Quantification

Learn more at: <https://www.biotechsupportgroup.com/category-s/335.htm>

NRicher™ Beads Are Versatile to A Variety of Bead Processing Formats

In addition to standard spin-filter formats, other formats compatible with the 50 µm NRicher™ beads are:

High Throughput Automation Compatible INTip™ SPE (DPX Technologies) Format

Aspirate and dispense
cycles mix NRicher™ beads
and solutions



The INTip™ SPE tip format improves ease of use and scalability to process multiple samples in parallel, utilizing 96-well plates and automated liquid handlers. The tip-based formats have been proven to be compatible with most automation platforms, i.e., Integra, Hamilton, etc. Please inquire for more information, as these formats are customized to the application and automation platform.

96-Well Vacuum or Pressure Filter Format

The NRicher™ beads can be readily processed in 96-well filter formats. Please inquire.



CONTACT US

We welcome your questions and comments regarding our products.

Tel: 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.

Email: sales@biotechsupportgroup.com

Web: www.biotechsupportgroup.com