



## BIOTECH SUPPORT GROUP

### SurfactAway™ SDS

#### SDS Detergent Removal

- Removes >99% detergent
- Very selective, virtually no cross-reactivity with other proteins
- Simple, just pipette, centrifuge and discard pellet
- Economical new surface technology, not based on hydrophobic chromatography

Detergents can often interfere with protein analysis. SurfactAway™ offers a simple and fast method to remove SDS and non-ionic detergents such as Triton. Recovery of protein is quantitative. **SurfactAway™ SDS** is especially designed for SDS removal and contains a precipitation buffer combined with a solid-phase binding suspension. **SurfactAway™ Non-ionic Triton**, is a sister product for non-ionic detergents. Both are applied in a simple protocol, just add, centrifuge and recover the protein solution.

### SDS-PAGE Electrophoresis



Lane 1: Plasma Control  
 Lane 2: Plasma Containing 0.1% SDS Treated With 0.25ml SurfactAway™ SDS  
 Lane 3: Plasma Containing 1.0% SDS Treated With 0.25ml SurfactAway™ SDS  
 Lane 4: IgG Control  
 Lane 5: IgG Containing 0.1% SDS Treated With 0.25ml SurfactAway™ SDS  
 Lane 6: IgG Containing 1.0% SDS Treated With 0.25ml SurfactAway™ SDS



## BIOTECH SUPPORT GROUP

Product	Size	# of Samples & Sample Size*	Item No.
<b>SurfactAway™ SDS</b>	30 ml	120, 1 ml samples	SA645-30
<b>SurfactAway™ SDS</b>	250 ml	1000, 1 ml samples	SA645-250
Product	Size	# of Samples & Sample Size*	Item No.
<b>SurfactAway™ Non-ionic</b>	30 ml	120, 1 ml samples	SA890-30
<b>SurfactAway™ Non-ionic</b>	250 ml	1000, 1 ml samples	SA890-250

### Storage

SurfactAway™ is supplied as an aqueous suspension of non-ionic adsorbent. Over time it will settle, and so before use, it should be shaken well to resuspend the solid-phase. Though it is stable at room temperature, suggested storage is at 4°C. Do not freeze. SurfactAway™ SDS should retain full activity when stored as directed for at least 6 months.

### Protocol

Actual free detergent concentration in biological samples can vary greatly, so the ratios shown are only intended to provide general guidance in use. As supplied, after centrifugation, the pellet occupies approximately 1/3 the suspension volume.

1. Resuspend SurfactAway™ SDS by gently shaking. Excessive shaking or mixing will cause foaming. It should be completely resuspended prior to use.
2. Add 1 ml of SurfactAway™ SDS to 4 ml of the sample. (1 : 4 ratio). Mix the sample by gently shaking periodically for 10 minutes.
3. Centrifuge sample at 16,000 G's for 1 minute - or - 1,000 G's for 15 minutes.
4. Decant supernatant containing macromolecules of interest and continue with purification.
5. Different sample volumes are easily scaled. Volume ratio can be adjusted up or down as required to remove the desired amount of detergent.

SurfactAway™ shows minimal cross reactivity with most serum proteins, but should not be used in excess.



## BIOTECH SUPPORT GROUP

### References

[Extraction and identification of electroimmunoprecipitated proteins from agarose gels.](#)

Journal of Immunological Methods Volume 330, Issues 1-2, 31 January 2008, Pages 24-33

### CONTACT US

**We welcome your questions and comments regarding our products.**

Tel: 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.

Fax: 732-274-2899

Email: [sales@biotechsupportgroup.com](mailto:sales@biotechsupportgroup.com)

Mail: 1 Deer Park Drive, Suite M, Monmouth Junction, NJ 08852

Web: [www.biotechsupportgroup.com](http://www.biotechsupportgroup.com)